**MTA ANGELUS**

# PORTUGUÊS

## INTRODUÇÃO

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## **ENGLISH**

### **INTRODUCTION**

MTA ANGELUS® is an bioceramic endodontic cement composed of several mineral oxides. It is constituted by thin hydrophilic particles. It is indicated in cases of root canal lateral and furcation perforations, internal resorption, reverse root filling, pulp capping, pulpotomy, apexification, and apexogenesis.

MTA (Mineral Trioxide Aggregate) provides the following benefits:

- Particle size which allows complete wetting during mixing (1, 2);
- Excellent marginal sealing; avoids penetration of tissue fluids in the root canal (3, 4, 5);
- Enclosing of root canal and furcation perforations through induction of periradicular cement formation (6, 2);
- Promotes the formation of a dentin bridge when used in pulp capping (3, 7);
- Unlike other cements, which demand a completely dry field, MTA is indicated when moisture control is inadequate (e.g., surgery for treatment of root perforation, reverse root filling), without loss of its properties (2).

### **COMPOSITION**

- Gray MTA: Tricalcium silicate, dicalcium silicate, tricalcium aluminate, calcium oxide, iron tetracalcium aluminate, bismuth oxide;
- White MTA: Tricalcium silicate, dicalcium silicate, tricalcium aluminate, calcium oxide, bismuth oxide.

### **PROPERTIES**

- Setting reaction: When mixed with distilled water it forms a gel that solidifies if kept in a wet environment. The initial setting time is approximately 10 minutes and the final is 15 minutes. It is not necessary to wait for the final set to continue treatment procedures (2);
- Hydrogen ion concentration (pH): After mixing, pH value is 10; in 3 hours it becomes highly alkaline (12)

(1, 8, 9);

- Radiopacity: Nearly matches that of gutta-percha. More radiopaque than dentine and bone (9);
- Compressive Strength: 40 MPa after 24 hours and 65 MPa after 21 days (9).

Sites of application do not receive direct occlusal load.

## **DIRECTIONS FOR USE**

1. Sterilize a glass slab, a metal spatula and all instruments for the insertion of MTA ANGELUS®;
2. Mix for 30 seconds the content of 1 sachet of MTA ANGELUS® (or 1 spoon of MTA ANGELUS®) with 1 drop of distilled water. The mixture should be homogeneous and with a consistency similar to wet sand;
3. Place the cement on the selected site with a sterilized amalgam carrier or other appropriate instrument;
4. Condense the cement with instruments such as amalgam condensers, a number 1 spatula or absorbent paper points moistened with distilled water.

**IMPORTANT:** If MTA ANGELUS® is not used immediately after mixing, its dehydration can be prevented and the working time increased by covering the mix on the glass slab with a wet gauze. If the mixed MTA dehydrates, it must be discarded.

## **INDICATIONS**

1. Treatment of perforations of root canal and furcation caused iatrogenically or by caries lesion (2) (Pics. 1 and 2);
2. Via canal treatment of root perforation due to internal resorption (2) (Pic. 3);
3. Surgical treatment of root perforation due to internal resorption (2) (Pic. 4);
4. Periapical surgery with reverse filling (5) (Pic. 5);
5. Pulp capping (7);
6. Pulpotomy (removal of affected coronal pulp to preserve vitality of remaining pulp tissue (2) (Pic. 6);
7. Apexogenesis (induction of root development in vital teeth with an inflamed coronal pulp) (2);
8. Apexification (induction of formation of a mineralized barrier at the root tip of young permanent teeth with incomplete root development and a necrotic pulp (2) (Pic. 7).

## **WARNINGS AND PRECAUTIONS**

- Use eye protection, mask and gloves when handling MTA ANGELUS®. In case of eye or skin contact, wash abundantly with water;
- Only open the sachet or bottle immediately prior to use. The MTA ANGELUS® powder is very sensitive to humidity. Close the bottle tightly after each use (2);
- Do not use MTA ANGELUS® to fill a root canal. Its viscosity is inadequate for this procedure and very difficult

to remove in case of a reentry;

- Do not use MTA ANGELUS® in areas of the tooth in contact with the gingival sulcus or it will be completely dissolved;
- Do not use MTA in areas with exposure to light in order to avoid darkening of the dental structure;
- Only use MTA ANGELUS® after remission of acute signs and symptoms of the endodontic disease. The acidic pH of endodontically compromised sites (lesions) prevents its setting reaction;
- Apply MTA ANGELUS® carefully. Similar to other endodontic cements, it is resorbed if extruded. However, the excess of any cement may impair the healing process;
- Do not use MTA on patients with reported sensitivity to any of its components.

**IMPORTANT:** The information provided in this manual is based on laboratory and clinical studies. The successful use of MTA ANGELUS® depends on a correct diagnosis, the operative technique, the condition of the treated tooth and the general health of the patient. This product must be used according to this manual.

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